In the Lab with Nuclear Scientists:

Critical and Subcritical Assemblies

American Nuclear Society Webinar Series on Educator Training

Jesse Norris Nuclear Criticality Safety Division Lawrence Livermore National Laboratory





Thermal Epithermal experiments (TEX)

- Critical experiments performed on the NCERC machines
 - Plutonium (²³⁹Pu and ²⁴⁰Pu)
 - Highly Enriched Uranium
 - Low-Temperature (-40°C)
- Based on a simple designs useful for modeling and validation
 - Design made to incorporate materials of interest
 - Tantalum, hafnium, lithium & more



Training Assembly for Criticality Safety (TACS)

 Subcritical assembly with Highly Enriched Uranium hemishells

Surrogate hemishells made of depleted uranium and steel

- Experiments for hands-on training include:
 - Approach to critical by
 - Fissile mass (with surrogates)
 - Reflector/moderator
 - Separation distance
 - Effect of reflection by hands
 - Effect of neutron poisons
 - Fissile material replacement measurements
- Operated on a hand-cranked
 Vertical Lift Machine



Inherently Safe Subcritical Assembly (ISSA)

- Subcritical assembly using MTR-type fuel assemblies
 - Aluminum surrogate assemblies
- Hands-on demonstrations
 - Opportunity for visitors to engage with fissile material
- Approach to critical by fissile mass and water height
 - Fixturing allows in-core diluents









Pulsed Neutron Die-Away (PNDA)

Nuclear experiments with or without fissile material

 Neutron generator pulsed to measure escaping neutrons

 Technique also applies to subcritical assemblies

Low-Scatter Facility

 Large room with false floor to reduce room return and background





Disclaimer

This document was prepared as an account of work sponsored by an agency of the United States government. Neither the United States government nor Lawrence Livermore National Security, LLC, nor any of their employees makes any warranty, expressed or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States government or Lawrence Livermore National Security, LLC. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States government or Lawrence Livermore National Security, LLC, and shall not be used for advertising or product endorsement purposes.